

## Portable Low-Cost Confocal Microscope Based on Pinhole Array and Smartphone

Presenter: Xuan-Ran, Wang(王炫然)

**Collaborator:**

Thesis Advisor: Prof. Ming-Che Chan 詹明哲

Thesis Advisor: Prof. Zu-Po Yang 楊斯博

ω *Photonic-Diagnosis LAB, NYCU, Taiwan*




# Introduction

## Advantage

- **Optical sectioning.**
- High resolution.
- Blocking defocused.
- Scattered light.
- Depth imaging.

## Disadvantage

- **Bulky system.**
- **High cost.**



US350,000=NTD10,000,000



# Conclusion

## limitation

The cost is approximately NTD 60,000 .

- Without ambient light
- Image need to be aligned
- Two smartphone


Components	NTD
Lens*4	2870
DMLP567	6000
FGL495	1300
EBS1	1150
Optomechanical components	22000
Smartphone*2	10600
Objective lens	14420
Total	58340



## Resolution

Horizontal Resolution of 20X is 2.19  $\mu\text{m}$ .  
Theoretical Resolution of 20X is 1.99  $\mu\text{m}$ .  
Vertical sensitivity of 20X is 0.5  $\mu\text{m}$ .




## Portable Confocal Microscope





# Confocal microscopy theory

## Pinhole array confocal microscope



# Master thesis Outline

- Introduction & Conclusion
- Theory
- Experimental process
- Results

# Experimental process


## Part1:Image process

Obtain a 2D image stack.

Align the transmission and reflection image stacks.

Select the region where the pinhole is clear, plot the z-axis curve to obtain the 1D longitudinal response curve.

Normalize each image stack, then subtract them to obtain the slope of the linear region.




## Part2:



Establish a calibration curve.

Divide the image's light intensity by the slope of the linear region to convert it into depth information.

Use MATLAB for 3D rendering to visualize the data and calculation results for analysis and interpretation.




Depth of the 2D image





# Experimental process



# Experimental process

Part2


Establish a calibration curve.

Divide the image's light intensity by the slope of the linear region to convert it into depth information.

Use MATLAB for 3D rendering to visualize the data and calculation results for analysis and interpretation.





Depth of 2D image:

$$\frac{\text{Differential intensity}}{\text{slope}(\frac{\text{Differential intensity}}{\text{depth}})}$$





-2.56E-2

MATLAB




Linear region



Slope:-2.56E-2  
(differential intensity/depth)




# Experimental results-Tunable lens



I(mA)	111.06	114.88	115.02	115.16	115.3	115.44	119.21
Focal length(mm)	76.49	75.125	75.07	75.02	74.978	74.93	73.65
Surface change(um)	-16.18um (2899.06)	-1.11um (2914.13)	-0.55um (2914.69)	0um (2915.24)	0.55um (2915.79)	1.1um (2916.34)	15.97um (2931.21)


(b) Curvature radius Slice 1

Slice2=74.93mm .....Slice6=75.125mm




(c)

Curvature radius



thick lens equation  
 $f = R/(n-1)$

Fitting



Residual

